Incomparable Dual Wave Light Quality

Dr's Light Clever

Model : CL-DC20

User manual

Eng

Kor

CE

Incomparable Dual Wave Light Quality

Dr's Light Clever

Model : CL-DC20

User manual

Eng

CE

1. Introduction	04
- Features	
- Application	
2. Safety information	05
- Symbols	
- Notes and reference	
- Installation and connection	
- Operation	
- Maintenance	
3. List of components	09
4. Parts and functions	10
5. Installation and connection	11
6. Direction for use	13
7. Selection of curing modes	14
8. Maintenance	15
- Cleaning and sterilization	
- Maintenance and replacement of the Battery	
9. Trouble shooting	16
10. Specification	17
11. Accessories	18
12. Warranty	18
13. Warranty card	19

1 Introduction

Thank you for purchasing the device.


The instructions described below have been designed to provide necessary information for successful operation, daily maintenance and safety for use of this unit.

Features

- Dual wavelength LED (peak 405nm & 460nm)
- Autoclaveable guide-cap (134°C)
- Variety of curing modes
- Easy-to read OLED screen
- Enhanced capacity battery
- Lithium-ion battery
- Interchangeable lithium-ion battery
- Chargeable battery alone
- Automatic sleep mode
- Ergonomic and familiar design

Application

- Polymerization of light-cure catalyst containing dental materials with a wavelength spectrum of 400 ~ 490nm.

2 Safety information

It is important to understand the following information to use this product safely.

Symbols

Indicates a potentially hazardous situation which, if not avoided, may result in serious injury or death.

Indicates a potentially hazardous situation which, if not avoided, may result in minor/moderate injury or in damage of property

Notes and References

Notes and References used in this manual are shown below.

Indicates that potential malfunctions may occur if the description is ignored.

Supplementary information to facilitate better operation.

2 Safety information

Installation and Connection

Warning

- The adaptor should be disconnected from the charger before connecting the power cable to prevent an electric shock.
- Keep away from sources of ignition to avoid fire or explosion.
- Protect the handpiece, charger, and battery from any liquid to avoid electric shock or functional error.
- Place the charger on a flat and stable surface. Do not drop handpiece or charger, in order to avoid internal defect that may cause overheating or fire.
- Do not place the electrical cord near a gas burner other sources of heat. Should the cord catch on fire, replace with a new one to avoid short-circuit.

Caution

- Do not use or leave the unit in areas of high temperature or exposed to direct sunlight; internal circuit defect may cause overheating or fire.
- Securely connect the dust-free power plug to the outlet. Use this device in well ventilated areas. Overheating may cause fire.
- Do not touch the terminal areas of charger or battery with fingers or conductive objects, such as metal wire or safety-pin. It may cause electric shock or fire.
- Do not use any charger other than supplied. It may result in electric shock or fire.

2 Safety information

Operation

Warning

- Do not use the device without guide-cap. It may result in burn cause by heat from LED.
- Do not plug the power cord with wet hands. This may result in an electric shock.
- In case of fluid leak from battery or partial discoloration/deformation of the Handpiece, immediately discontinue use.
- Place the battery in correct direction when charging. Wrong directed placement may cause malfunction, explosion, or fire.
- In case of partial discoloration of the charger, immediately disconnect the power plug and repair the charger.
- In case battery fluid contacts the eye, immediately flush with copious amount of clean water and seek medical attention.
- Do not use this product on patients with a history of cataract, retinal damage or allergic reaction to light.
- Do not use this product on a person with a pacemaker as it may cause malfunction of the pacemaker.
- Avoid exposing the light on intraoral soft tissue as it may affect the soft tissue.
- Use carefully the guide-cap to avoid broken lens on the guide-cap.
- Do not put the battery on the charger upside down. This may cause an explosion or fire.

Caution

- Disconnect the power cord by grasping the power plug to avoid injuries or burns.
- Avoid damage or use of unnecessary force on the power cord. This may cause an electric shock or fire.
- Should battery fluid contact the skin or clothes, immediately wash it off with copious amounts of water to avoid skin damage.
- Take the battery out of the handpiece if it won't be used for more than 1 week.
- Connect the power plug only when battery is being charged. Insulation degradation may cause electric shock or fire.
- Do not use this product for any purposes other than dental use.

2 Safety information

Maintenance

Warning

- Use only the battery supplied to avoid fluid leak or explosion.
- Never disassemble or overhaul the product. Unusual performance may cause injuries, electric shock, or fire.
- Ensure that charger and battery are always separated from the handpiece prior to cleaning.


Caution


- If the OLED screen malfunctions, remove the battery from the handpiece and reinsert it.
- User must be in charge of operation, maintenance and check.
- Take care a lens on the guide-cap to prevent broken when cleaning.

3 List of components


Please check the contents described below when removing the contents from the original packaging.


Handpiece


Charger


Battery


Guide-cap (x 3)

Adapter


Power cable

User manual

4 Parts and functions

① Handpiece

- ① Start button : Operation and turn On/Off
- ② Mode button : Curing mode select
- ③ Control butto : Curing time select
- ④ OLED screen : Display device state

② Guide-cap: Autoclaveable, easy to replace

③ Battery

④ Charger: Recharge battery display charging process

⑤ Adapter


⑥ Power cable

5 Installation and connection

Refer to the information supplied on safety precautions during installation and connection.

1. Connection between adapter and charger

- 1) Securely insert the power cable with 100~240V power outlet.
- 2) Connect power cable with adapter.
- 3) Connect adapter with charger


Note

- Use only the original adapter supplied.
- When the adapter is connected with charger, red and blue light will blink.
- If red and blue light is not blinking when adapter is connected, it is a charger malfunction.

2. Inserting battery


- ① Carefully insert the battery into the handpiece as directed by the arrow.
- ② The battery turns with the arrow, as seen in image.
***Avoid forcing the battery.**

5 Installation and connection

3. Inserting guide-cap

- Slowly insert the guide-cap into the handpiece as picture. Insert it till right position.


Reference

- Replace the guide-cap when it is dirty, or intensity is lower than normal condition caused by scratch or broken on lens.
- For convenience, sterilize the guide-cap before use.

4. Charging

- Insert the handpiece into the vertical slot at the front of the charger.
- Red light on : Charging in progress
- Blue light on : Charging process is complete

Note

- If the remaining battery is extremely low, red light may flicker while recharging, which is not a malfunction.


Reference

- It will take approximately 3~4 hours before the battery will be fully charged.


6 Direction for use

1. Activating handpiece

1) Press any button, then OLED screen will be on standby mode and time display used last.

2) Select curing mode and time setting with mode button and control button.

Initial screen

Standby screen
(last used mode & time)

Reference

- When battery inserted, the unit is automatically turned on.

2. Start of exposure

Press start button to start exposure with the selected mode.

Reference

- Press any button if you wish to end the curing process.
- Alarm during operating informs time passage.

3. Battery icon

This icon indicates the power remaining in the battery. Recharge or replace the battery when "LOW BATTERY" sign appears on the display screen.

Reference


- Sleep mode : If unit is left unused for a few minutes, all functions are automatically transferred to sleep mode with a short beep. Hold down start button for approximately 3-5 seconds to shut off light. If you want to turn on, press any button.
- Set-up memory : The handpiece will record the setting used last.

7 Selection of curing mode

Curing mode can be set by pressing the mode button.
Various curing programs are available in this unit.


1. High mode

Time setting: 5, 10 sec.
Quick polymerizing with high intensity


2. Low mode

Time setting: 10, 20, 30 sec.
Polymerizing with low intensity


3. Soft start mode

Time setting: 20 sec. (fixed)
Special mode to avoid polymerization shrinkage

4. Pulse mode

Time setting: 10, 20 sec.
Special mode to avoid polymerization shrinkage

Cleaning and sterilization

- Clean the handpiece and charger gently with clean, soft, lint-free cloth.
- Avoid using solvents which contain ethanol or acetone.
- Do not use sharp or pointed tools to clean as they may damage or scratch the plastic.
- Clean the guide-cap with a soft cloth and disinfect. The guide-cap is autoclaveable.
- Remove restorative materials stuck to the guide-cap gently with a plastic spatula. Avoid using sharp tool as they may scratch the surface of the lens and affect the intensity of light.

Maintenance and replacement of the battery

- Exposure to high room temperature may shorten the battery life and affect device performance. Use the device in normal temperature.
- It is recommended for rechargeable battery to be charged immediately before use to maximize efficiency.
- Remove the battery from the handpiece, if not being used for a long time.
- If the running time declines with a full-charged battery, replace the battery before use.

9 Trouble shooting

If your curing light unit is not operating normally, please refer to the following points.

	Trouble	Cause and solution
1	Control button does not work	If the control button is pressed with excess force, it will not work. In such case press the button gently.
2	Display screen does not work	If the battery is completely discharged or connected incorrectly, display will not appear. Recharge the Battery or reconnect the Battery to reset the device.
3	Power off while in operation	If the start button is pressed for a moment, power will be turned off. Press any button again. Check the battery connection again. If the battery is completely discharged, device can be turned off. Recharge the battery.
4	Not able to charge	If the power cord or etc. is not correctly connected, the battery can not be charged. If the battery or handpiece is not properly placed on the charger, the battery can not be charged.
5	Poor intensity	Dirty, damage of guide-cap can make poor intensity. Clean the guide-cap or replace new one.

* If the trouble is not settled after upper solutions, please contact local dealer or Good Doctors Co., Ltd.

10 Specification

Handpiece	Operation voltage	3.7V
	Peak wavelength	405nm & 460nm
	Dimensions (included battery)	Approx. 232 x 40 x 26mm
	Weight (included battery)	Approx. 130g

Charger	Operation voltage	5V
	Dimensions	Approx. 110 x 105 x 68mm
	Weight	Approx. 100g


Battery	Type	Lithium Ion battery
	Nominal voltage	DC 3.7V
	Dimensions	Approx. 75 x 25 x 25mm
	Weight	Approx. 52g

Guide-cap	Dimensions	Approx. 64 x 16 x 16mm
	Weight	Approx. 5g


Adapter	Input voltage	AC 100 ~ 240V
	Output voltage	DC 5V

* This specification is subject to change without prior notice

11 Accessories

Guide-cap (x 3)

Battery

12 Warranty

This product is shipped after strict inspection. Should the product have any problem under normal usage within one year from the date of purchase, repair will be free of charge.

However, repair of consumables (guide-cap, battery, etc.) are not covered by the warranty. For more information, refer to the warranty card attached.

EC REP Gruneworld GmbH
Georg-Weiss-Weg 9, 53125 Bonn, Germany
TEL: +49(0)288/7077695 FAX: +49(0)288/7077696
Website: www.gooddrs.de

 Good Doctors Co., Ltd.
#208, B-Dong, 283 Bupyeong-daero, Bupyeong-gu,
Incheon, Korea (Woolim Lions Valley, Cheongcheon-dong)
(Zip code : 403-911)
Tel : + 82 32 424 6325 Fax : + 82 32 424 6326
E-mail : info@gooddrs.com Website : www.gooddrs.com

Good Doctors Good Doctors Good Doctors Good Doctors Good Doctors

We guarantee the quality of the product as follows

< Assurance details >

1. If there is any defect of manufacturer or product has any fault itself within the assurance period, it will be replaced with the same product.
 - Malfunction and damage according to user's carelessness,
 - Malfunction and damage according to natural disaster such as fire, earthquake and flood etc.
 - Malfunction and damage according to user's intention or accident.
 - Malfunction and damage according to using incompatible product.
2. Below cases are not effective for the free refund.
 - Malfunction and damage according to using incompatible product.
3. The manufacturer does not accept any liability resulting misuse.

Name of product :

Name of model :

Customer information

- Office name :
- Office phone :

Distributor

- Dealer name :

Assurance period

- Date of purchase : from _____ , 20____
- Product is for 1 year except accessories.

Good Doctors Co., Ltd.
 #208, B-Dong, 283 Bupyeong-daero, Bupyeong-gu, Incheon, Korea (Woolim Lions Valley, Cheongcheon-dong) (Zip code : 403-911)
 Tel: +82 32 424 6325 Fax: +82 32 424 6326

Good Doctors Good Doctors Good Doctors Good Doctors Good Doctors

Incomparable Dual Wave Light Quality

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