clinic line

ventura impress 2

Addition putty silicone of soft consistency

Features & benefits:

High resistance to deformation.

High tear strength.

Excellent elastic recovery.

Specially indicated for:

One and two step putty-wash technique for crown & bridge work inlays and onlays impressions.



Addition putty silicone of soft consistency. Avalable in regular and fast setting times.



Available in:

Refills:

2003 - ventura impress 2 putty soft regular set - 300 + 300 ml. 2016 - ventura impress 2 putty soft fast set - 300 + 300 ml.

Properties:

	Putty soft	Putty sof	
	regular Set	fast set	
Working time:	2'.	1'15".	
Time in the oral cavit	y: 2'30".	1'30".	
Setting time:	3'.	2'.	
Shore A hardness:	60.	60.	
Colours:	Yellow.	Blue.	

