

ventura impress 2

light body

Light consistency addition silicone

Features & benefits:

- Precision in detail reproduction.
- Thixotropic.
- High resistance to deformation.
- Very high tear strength.
- Excellent elastic recovery.
- Extremely hydrophilic.



Light consistency addition silicone for accurate impressions.

Specially indicated for:

One and two step putty-wash technique for crown & bridge work inlays and onlays impressions. Also for reline impressions.



Available in:

Refills:

- 2004 - ventura impress 2 light body fast set - 2 x 50 ml.
- 2018 - ventura impress 2 light body fast set - 4 x 50 ml.
- 2017 - ventura impress 2 light body regular set - 2 x 50 ml.
- 2019 - ventura impress 2 light body regular set - 4 x 50 ml.

Accessories:

- 2010 - ventura mixing tips yellow - 50 u.
- 2011 - ventura intraoral tips yellow - 50 u.

Properties:

	Light body regular set	Light body fast set
Working time :	2'15"	1'30"
Time in the oral cavity :	3'	1'30"
Setting time:	3'30"	2'
Shore A hardness:	53.	48.
Colours:	Green.	Orange.

